

LSC as Powerful and Fast Tool for In-situ Measurement of Natural Radionuclides in Water

T. L. Möbius*, E. Frenzel**, J. Haaslahti°, K. Kamolchote°°

* Institut National des Sciences et Techniques Nucléaires, INSTN, Antananarivo, Madagascar

** FCI Frenzel Consulting & Instruments, Straubenhardt, Germany

° HIDEX, Turku, Finland

°° Silpakorn University, Nakhon Pathom, Thailand

Revised version of a presentation, given at Conference BAER in Gaithersburg, USA, October 1999

Acknowledgement

The authors are grateful to Dr.S. Moebius, Research Centre of Karlsruhe, Germany. Data shown in this paper are results from comprehensive investigations coordinated by Siegurd Möbius, Research Center Karlsruhe. They have been presented at the Conference BAER in Gaithersburg, October 1999 and will be published in 'Radioactivity & Radiochemistry'.

Abstract

When natural radionuclides have to be determined in water samples sampling and sample conservation during transport to the laboratory generally is considered as main cause for errors. Especially when Radon and its Progenies are investigated on site analysis is advantageous. Various measuring principles using different nuclides for the estimation of the Radon level aggravate the assessment of results.

The paper describes methods for in-situ measurement and enrichment of Radium, Radon and daughter nuclides followed by in-situ Liquid Scintillation Spectrometry with the portable liquid scintillation system Triathler. The methods for Radon include extraction into an organic cocktail like Beta Plate Scint or those based on toluene as solvent and gel counting. Radium is enriched by adsorption on MnO₂ coated polyamide granulates or on a BaSO₄ loaded cationic resin. These procedures have shown to provide reliable results through various intercomparison runs (e.g. Thale in Germany 1999, Luxembourg 1997). Data are presented and compared with other techniques.

When both a- and b-spectra are evaluated a survey on the state of radioactive equilibrium is available. Using the a-channel from a-/b-pulse shape discrimination a lower limit of detection of 0.05 Becquerel Radon per Liter has been calculated for a 100 min counting period with the Triathler™. In addition a two-dimensional spectrometric plot of pulse length versus pulse height provides a powerful quality control of counting data. Preliminary results from field measurements are discussed.

Product Information

Finland and International Hidex Oy tel +358-2-2750557 fax +358-2-2410075 info@hidex.com **Canada** Gamble Technologies Limited Tel: 905 812 9200 info@gtl.ca **France** ScienceTec Tel +3316446 2400 info@sciencetec.fr **Germany (Austria, Czech Republic)** FCI Frenzel Consulting & Instruments Tel +49 7082 948 290 eberhardfrenzel@aol.com **Hong Kong** Line Analytics Life Sciences Ltd. Tel + 852 2578 5839 line@lineanalytics.com **Japan** Sowa Trading Co Inc. Tel +81 3 3862 2700 tokyo@sowa-trading.co.jp **Korea** SIMS Co.Ltd Tel +82 2 576 9727 dmlee@sims.com **Russia** Pribori Oy c/o Wallac Oy Tel +358 2 267 8111 maria.kivisto@pribori.pp.fi **Scandinavia:** Reaction-Lab Oy Tel +358 2 410 1746 info@reaction-lab.com **Spain** Izasa s.a. Tel +34 91 6630500 jbruna@izasa.es **Switzerland** Uehlinger Pfiffner AG Tel +41 62 7212281 mail@uehlinger-pfiffner.ch **United Kingdom** Southern Scientific Ltd. Tel +44 1903 6044000 info@southern-scientific.co.uk **U.S.A** Bioscan Inc Tel + 1 202 338 0974 sales@bioscan.com

www.hidex.com